



Physician Name and Facility

Physician Name			Patient Name and DOB		
Date Ordered			Case Accession Number		
☐	Prepare H&E ONLY (Technical)	☐	Surgical Consultation		
☐	Gross & Prepare H&E (Technical)	☐	Surgical Sign Out ONLY (Pro Fee)		
☐	Prepare H&E w/Interp (Pro and Tech)	☐	Notification of pending Flow to ARL		
☐	IHC/Specialty Stain (Technical ONLY)	☐	IHC/Specialty Stain (Pro and Tech)		

Epithelial		Hematologic		Other		Neuroendocrine		Predictive Markers	
☐	Calrentinin	☐	BCL-2	☐	AFP	☐	Calcitonin	☐	ER (ErR)
☐	CEA/Polyclonal	☐	BCL-6	☐	CDX-2	☐	CD56	☐	PR (PgR)
☐	EMA	☐	CD3	☐	DOG1	☐	Chromogranin A	☐	Her-2/Neu
☐	Pan Keratin (AE1/AE3)	☐	CD5	☐	Fascin	☐	Synaptophysin	☐	MMR
☐	OSCAR (Cytokeratin)	☐	CD10	☐	GATA3	☐	Thyroglobulin		
☐	CK5/6	☐	CD15	☐	GFAP				
☐	CK7	☐	CD20 (L26)	☐	HSA	Soft Tissue		Recut/Deeper	
☐	CK20	☐	CD23	☐	H.pylori	☐	Actin (SM)	☐	Recut(s)
☐	p16	☐	CD30	☐	Mammaglobin	☐	CD34	☐	Deeper(s)
☐	p40	☐	CD45 (LCA)	☐	Ki-67	☐	Desmin	☐	Cut through block
☐	E-cadherin	☐	CD68	☐	p53	☐	Factor VIII		
☐	EP-CAM (Ber4)	☐	CD117	☐	PLAP	☐	STAT6		
☐	PAX8	☐	CD138	☐	PSA				
☐	Napsin A	☐	Cyclin D1	☐	SALL4	Melanocytic		Quality Assurance:	
☐	NKX3.1	☐	Kappa	☐	Tryptase	☐	HMB-45	HT/HTL Initials: _____	
		☐	Lambda	☐	TTF-1	☐	MART-1	Pathologist Initials: _____	
SS Group I		☐	MPO	☐	Vimentin	☐	S100	Controls: Satisfactory, unless otherwise stated. Satisfactory means that positive control(s) stained appropriately as well as patient tissue appears technically acceptable.	
☐	AFB	☐	MUM1	☐	CMV	☐	SOX10		
☐	B&B Gram	☐	PAX5						
☐	GMS	☐	C-MYC						
☐	PAS/Fungus				IHC Panels				
					☐	NSCLC (p40, TTF1)			
					☐	Hodkin's (CD3, 15, 20, 30, 45, PAX5, Fascin, MUM1)			
SS Group II		☐	PAS	☐	Germ Cell Tumor (AFP, CD30, CD117, OSCAR)				
☐	Alcian Blue pH2.5	☐	PAS w/Diastase	☐	Diffuse Large B-Cell (CD3, 5, 10, 20, 23, 30, BCL2, Bcl6, CyclinD1, MUM1, Ki-67, cMYC)				
☐	Alkaline Congo Red	☐	G&S Reticulin	☐	Mesothelioma (AE1/AE3, EMA, Desmin, Calrentin, CK5/6, TTF1, CEA(P), EP-CAM(Ber4))				
☐	Gomori's Fe	☐	Masson Trichrome						
☐	Mayer's Mucicarmine	☐	VVG-Elastin	IHC and Stain Disclaimer: Use of fixatives(s) other than 10% NBF may not yield equivalent					

Controls: Satisfactory, unless otherwise stated.
Satisfactory means that positive control(s) stained appropriately as well as patient tissue appears technically acceptable.

IHC and Stain Disclaimer: Use of fixatives(s) other than 10% NBF may not yield equivalent or satisfactory results. Our IHC and specialty stains ONLY have been validated on 10% NBF fixed, non-decalcified tissues.